

Antibodies Product Datasheet

Product Name:	MC2-R Antibody
Catalog Number:	G210
Reactivity:	Human, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ACTHR. AA range:248-297
Concentration:	1 mg/ml
Gene Name:	MC2R
Protein Name:	Adrenocorticotrophic hormone receptor
Other Name:	MC2R; ACTHR; Adrenocorticotrophic hormone receptor; ACTH receptor; ACTH-R; Adrenocorticotropin receptor; Melanocortin receptor 2; MC2-R
Human Swiss Prot No.:	Q01718
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MC2-R Polyclonal Antibody detects endogenous levels of MC2-R protein.
Storage and Stability:	-20°C/1 year