

Antibodies Product Datasheet

Product Name:	AT1 Antibody
Catalog Number:	G212
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AGTR1. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	AGTR1
Protein Name:	Type-1 angiotensin II receptor
Other Name:	AGTR1; AGTR1A; AGTR1B; AT2R1; AT2R1B; Type-1 angiotensin II receptor; AT1AR; AT1BR; Angiotensin II type-1 receptor; AT1
Human Swiss Prot No.:	P30556
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AT1 Polyclonal Antibody detects endogenous levels of AT1 protein.
Storage and Stability:	-20°C/1 year