

Antibodies Product Datasheet

Product Name:	BAI-3 Antibody
Catalog Number:	G215
Reactivity:	Human, Mouse
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BAI3. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	BAI3
Protein Name:	Brain-specific angiogenesis inhibitor 3
Other Name:	BAI3; KIAA0550; Brain-specific angiogenesis inhibitor 3
Human Swiss Prot No.:	O60242
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BAI-3 Polyclonal Antibody detects endogenous levels of BAI-3 protein.
Storage and Stability:	-20°C/1 year