

Antibodies Product Datasheet

Product Name:	Bradykinin B1 R Antibody
Catalog Number:	G216
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human BDKRB1. AA range:201-250
Concentration:	1 mg/ml
Gene Name:	BDKRB1
Protein Name:	B1 bradykinin receptor
Other Name:	BDKRB1; BRADYB1; B1 bradykinin receptor; B1R; BK-1 receptor
Human Swiss Prot No.:	P46663
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Bradykinin B1 R Polyclonal Antibody detects endogenous levels of Bradykinin B1 R protein.
Storage and Stability:	-20°C/1 year