

Antibodies Product Datasheet

Product Name:	CGI-43 Antibody
Catalog Number:	G223
Reactivity:	Human, Mouse, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CG028. AA range:411-460
Concentration:	1 mg/ml
Gene Name:	CCZ1/C7orf28B
Protein Name:	Vacuolar fusion protein CCZ1 homolog/Vacuolar fusion protein CCZ1 homolog-like
Other Name:	CCZ1; C7orf28A; CGI-43; Vacuolar fusion protein CCZ1 homolog; CCZ1B; C7orf28B; Vacuolar fusion protein CCZ1 homolog B; Vacuolar fusion protein CCZ1 homolog-like
Human Swiss Prot No.:	P86791/P86790
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CGI-43 Polyclonal Antibody detects endogenous levels of CGI-43 protein.
Storage and Stability:	-20°C/1 year