

Antibodies Product Datasheet

Product Name:	CysLTR1 Antibody
Catalog Number:	G224
Reactivity:	Human, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CLTR1. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	CYSLTR1
Protein Name:	Cysteinyl leukotriene receptor 1
Other Name:	CYSLTR1; CYSLT1; Cysteinyl leukotriene receptor 1; CysLTR1; Cysteinyl leukotriene D4 receptor; LTD4 receptor; G-protein coupled receptor HG55; HMTMF81
Human Swiss Prot No.:	Q9Y271
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CysLTR1 Polyclonal Antibody detects endogenous levels of CysLTR1 protein.
Storage and Stability:	-20°C/1 year