

Antibodies Product Datasheet

Product Name:	CB2 Antibody
Catalog Number:	G227
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CNR2. AA range:191-240
Concentration:	1 mg/ml
Gene Name:	CNR2
Protein Name:	Cannabinoid receptor 2
Other Name:	CNR2; Cannabinoid receptor 2; CB-2; CB2; hCB2; CX5
Human Swiss Prot No.:	P34972
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CB2 Polyclonal Antibody detects endogenous levels of CB2 protein.
Storage and Stability:	-20°C/1 year