

Antibodies Product Datasheet

Product Name:	EDG-3 Antibody
Catalog Number:	G235
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EDG3. AA range:115-164
Concentration:	1 mg/ml
Gene Name:	S1PR3
Protein Name:	Sphingosine 1-phosphate receptor 3
Other Name:	S1PR3; EDG3; Sphingosine 1-phosphate receptor 3; S1P receptor 3; S1P3; Endothelial differentiation G-protein coupled receptor 3; Sphingosine 1-phosphate receptor Edg-3; S1P receptor Edg-3
Human Swiss Prot No.:	Q99500
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EDG-3 Polyclonal Antibody detects endogenous levels of EDG-3 protein.
Storage and Stability:	-20°C/1 year