

Antibodies Product Datasheet

Product Name:	EDG-4 Antibody
Catalog Number:	G236
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EDG4. AA range:271-320
Concentration:	1 mg/ml
Gene Name:	LPAR2
Protein Name:	Lysophosphatidic acid receptor 2
Other Name:	LPAR2; EDG4; LPA2; Lysophosphatidic acid receptor 2; LPA receptor 2; LPA-2; Lysophosphatidic acid receptor Edg-4
Human Swiss Prot No.:	Q9HBW0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EDG-4 Polyclonal Antibody detects endogenous levels of EDG-4 protein.
Storage and Stability:	-20°C/1 year