

Antibodies Product Datasheet

Product Name:	EDG-6 Antibody
Catalog Number:	G238
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EDG6. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	S1PR4
Protein Name:	Sphingosine 1-phosphate receptor 4
Other Name:	S1PR4; EDG6; Sphingosine 1-phosphate receptor 4; S1P receptor 4; S1P4; Endothelial differentiation G-protein coupled receptor 6; Sphingosine 1-phosphate receptor Edg-6; S1P receptor Edg-6
Human Swiss Prot No.:	O95977
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EDG-6 Polyclonal Antibody detects endogenous levels of EDG-6 protein.
Storage and Stability:	-20°C/1 year