

Antibodies Product Datasheet

Product Name:	EDG-7 Antibody
Catalog Number:	G239
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EDG7. AA range:281-330
Concentration:	1 mg/ml
Gene Name:	LPAR3
Protein Name:	Lysophosphatidic acid receptor 3
Other Name:	LPAR3; EDG7; LPA3; Lysophosphatidic acid receptor 3; LPA receptor 3; LPA-3; Lysophosphatidic acid receptor Edg-7
Human Swiss Prot No.:	Q9UBY5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	EDG-7 Polyclonal Antibody detects endogenous levels of EDG-7 protein.
Storage and Stability:	-20°C/1 year