

Antibodies Product Datasheet

Product Name:	ETAR Antibody
Catalog Number:	G241
Reactivity:	Human, Mouse, Monkey
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human EDNRA. AA range:378-427
Concentration:	1 mg/ml
Gene Name:	EDNRA
Protein Name:	Endothelin-1 receptor
Other Name:	EDNRA; ETA; ETRA; Endothelin-1 receptor; Endothelin A receptor; ET-A; ETA-R; hET-AR
Human Swiss Prot No.:	P25101
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ETAR Polyclonal Antibody detects endogenous levels of ETAR protein.
Storage and Stability:	-20°C/1 year