

Antibodies Product Datasheet

Product Name:	GPR37L1 Antibody
Catalog Number:	G247
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human ETBR2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	GPR37L1
Protein Name:	Endothelin B receptor-like protein 2
Other Name:	GPR37L1; ETBRLP2; Endothelin B receptor-like protein 2; ETBR-LP-2; G-protein coupled receptor 37-like 1
Human Swiss Prot No.:	O60883
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR37L1 Polyclonal Antibody detects endogenous levels of GPR37L1 protein.
Storage and Stability:	-20°C/1 year