

Antibodies Product Datasheet

Product Name:	FPRL2 Antibody
Catalog Number:	G252
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human FPRL2. AA range:304-353
Concentration:	1 mg/ml
Gene Name:	FPR3
Protein Name:	N-formyl peptide receptor 3
Other Name:	FPR3; FPRH1; FPRL2; N-formyl peptide receptor 3; FMLP-related receptor II; FMLP-R-II; Formyl peptide receptor-like 2
Human Swiss Prot No.:	P25089
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FPRL2 Polyclonal Antibody detects endogenous levels of FPRL2 protein.
Storage and Stability:	-20°C/1 year