

Antibodies Product Datasheet

Product Name:	GPR137C Antibody
Catalog Number:	G264
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR137C. AA range:256-305
Concentration:	1 mg/ml
Gene Name:	GPR137C
Protein Name:	Integral membrane protein GPR137C
Other Name:	GPR137C; TM7SF1L2; Integral membrane protein GPR137C; Transmembrane 7 superfamily member 1-like 2 protein
Human Swiss Prot No.:	Q8N3F9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR137C Polyclonal Antibody detects endogenous levels of GPR137C protein.
Storage and Stability:	-20°C/1 year