

Antibodies Product Datasheet

Product Name:	GALR3 Antibody
Catalog Number:	G269
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GALR3. AA range:291-340
Concentration:	1 mg/ml
Gene Name:	GALR3
Protein Name:	Galanin receptor type 3
Other Name:	GALR3; GALNR3; Galanin receptor type 3; GAL3-R; GALR-3
Human Swiss Prot No.:	O60755
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GALR3 Polyclonal Antibody detects endogenous levels of GALR3 protein.
Storage and Stability:	-20°C/1 year