

Antibodies Product Datasheet

Product Name:	GHRH-R Antibody
Catalog Number:	G271
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GHRHR. AA range:351-400
Concentration:	1 mg/ml
Gene Name:	GHRHR
Protein Name:	Growth hormone-releasing hormone receptor
Other Name:	GHRHR; Growth hormone-releasing hormone receptor; GHRH receptor; Growth hormone-releasing factor receptor; GRF receptor; GRFR
Human Swiss Prot No.:	Q02643
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GHRH-R Polyclonal Antibody detects endogenous levels of GHRH-R protein.
Storage and Stability:	-20°C/1 year