

Antibodies Product Datasheet

Product Name:	GPR126 Antibody
Catalog Number:	G290
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR126. AA range:431-480
Concentration:	1 mg/ml
Gene Name:	GPR126
Protein Name:	G-protein coupled receptor 126
Other Name:	GPR126; DREG; VIGR; G-protein coupled receptor 126; Developmentally regulated G-protein-coupled receptor; Vascular inducible G protein-coupled receptor
Human Swiss Prot No.:	Q86SQ4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR126 Polyclonal Antibody detects endogenous levels of GPR126 protein.
Storage and Stability:	-20°C/1 year