

Antibodies Product Datasheet

Product Name:	GPRC5B Antibody
Catalog Number:	G318
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPRC5B. AA range:61-110
Concentration:	1 mg/ml
Gene Name:	GPRC5B
Protein Name:	G-protein coupled receptor family C group 5 member B
Other Name:	GPRC5B; RAIG2; G-protein coupled receptor family C group 5 member B; A-69G12.1; Retinoic acid-induced gene 2 protein; RAIG-2
Human Swiss Prot No.:	Q9NZH0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPRC5B Polyclonal Antibody detects endogenous levels of GPRC5B protein.
Storage and Stability:	-20°C/1 year