

Antibodies Product Datasheet

Product Name:	GPRC5C Antibody
Catalog Number:	G319
Reactivity:	Human, Mouse, Rat, Monkey
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPRC5C. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	GPRC5C
Protein Name:	G-protein coupled receptor family C group 5 member C
Other Name:	GPRC5C; RAIG3; PSEC0087; G-protein coupled receptor family C group 5 member C; Retinoic acid-induced gene 3 protein; RAIG-3
Human Swiss Prot No.:	Q9NQ84
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPRC5C Polyclonal Antibody detects endogenous levels of GPRC5C protein.
Storage and Stability:	-20°C/1 year