

Antibodies Product Datasheet

Product Name:	GPR31 Antibody
Catalog Number:	G334
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR31. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	GPR31
Protein Name:	12-(S)-hydroxy-5,8,10,14-eicosatetraenoic acid receptor
Other Name:	GPR31; 12-(S)-hydroxy-5; 8,10,14-eicosatetraenoic acid receptor; 12-(S)-HETE receptor; 12-HETER; G-protein coupled receptor 31
Human Swiss Prot No.:	O00270
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR31 Polyclonal Antibody detects endogenous levels of GPR31 protein.
Storage and Stability:	-20°C/1 year