

Antibodies Product Datasheet

Product Name:	GPR37 Antibody
Catalog Number:	G339
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR37. AA range:211-260
Concentration:	1 mg/ml
Gene Name:	GPR37
Protein Name:	Probable G-protein coupled receptor 37
Other Name:	GPR37; Probable G-protein coupled receptor 37; Endothelin B receptor-like protein 1; ETBR-LP-1; Parkin-associated endothelin receptor-like receptor; PAELR
Human Swiss Prot No.:	O15354
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR37 Polyclonal Antibody detects endogenous levels of GPR37 protein.
Storage and Stability:	-20°C/1 year