

Antibodies Product Datasheet

Product Name:	GPR62 Antibody
Catalog Number:	G348
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR62. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	GPR62
Protein Name:	Probable G-protein coupled receptor 62
Other Name:	GPR62; Probable G-protein coupled receptor 62; G-protein coupled receptor GPCR8; hGPCR8; G-protein coupled receptor KPG_005
Human Swiss Prot No.:	Q9BZJ7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR62 Polyclonal Antibody detects endogenous levels of GPR62 protein.
Storage and Stability:	-20°C/1 year