

Antibodies Product Datasheet

Product Name:	GPR63 Antibody
Catalog Number:	G349
Reactivity:	Human, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR63. AA range:370-419
Concentration:	1 mg/ml
Gene Name:	GPR63
Protein Name:	Probable G-protein coupled receptor 63
Other Name:	GPR63; PSP24B; Probable G-protein coupled receptor 63; PSP24-2; PSP24-beta
Human Swiss Prot No.:	Q9BZJ6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR63 Polyclonal Antibody detects endogenous levels of GPR63 protein.
Storage and Stability:	-20°C/1 year