

Antibodies Product Datasheet

Product Name:	GPR82 Antibody
Catalog Number:	G354
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GPR82. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	GPR82
Protein Name:	Probable G-protein coupled receptor 82
Other Name:	GPR82; Probable G-protein coupled receptor 82
Human Swiss Prot No.:	Q96P67
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR82 Polyclonal Antibody detects endogenous levels of GPR82 protein.
Storage and Stability:	-20°C/1 year