

Antibodies Product Datasheet

Product Name:	mGluR1 Antibody
Catalog Number:	G363
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	Immunohistochemistry: 1/100 - 1/300. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human GRM1. AA range:251-300
Concentration:	1 mg/ml
Gene Name:	GRM1
Protein Name:	Metabotropic glutamate receptor 1
Other Name:	GRM1; GPRC1A; MGLUR1; Metabotropic glutamate receptor 1; mGluR1
Human Swiss Prot No.:	Q13255
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	mGluR1 Polyclonal Antibody detects endogenous levels of mGluR1 protein.
Storage and Stability:	-20°C/1 year