

Antibodies Product Datasheet

Product Name:	Latrophilin-1 Antibody
Catalog Number:	G377
Reactivity:	Human, Mouse, Rat
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LPHN1. AA range:561-610
Concentration:	1 mg/ml
Gene Name:	LPHN1
Protein Name:	Latrophilin-1
Other Name:	LPHN1; KIAA0821; LEC2; Latrophilin-1; Calcium-independent alpha-latrotoxin receptor 1; CIRL-1; Lectomedin-2
Human Swiss Prot No.:	O94910
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Latrophilin-1 Polyclonal Antibody detects endogenous levels of Latrophilin-1 protein.
Storage and Stability:	-20°C/1 year