

Antibodies Product Datasheet

Product Name:	Latrophilin-2 Antibody
Catalog Number:	G378
Reactivity:	Human, Mouse, Rat
Applications:	IHC-p, IF/ICC, ELISA
Dilutions:	WB 1:500-2000 ,Immunohistochemistry: 1/100 - 1/300. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LPHN2. AA range:551-600
Concentration:	1 mg/ml
Gene Name:	LPHN2
Protein Name:	Latrophilin-2
Other Name:	LPHN2; KIAA0786; LEC1; LPHH1; Latrophilin-2; Calcium-independent alpha-latrotoxin receptor 2; CIRL-2; Latrophilin homolog 1; Lectomedin-1
Human Swiss Prot No.:	O95490
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Latrophilin-2 Polyclonal Antibody detects endogenous levels of Latrophilin-2 protein.
Storage and Stability:	-20°C/1 year