

Antibodies Product Datasheet

Product Name:	LHR Antibody
Catalog Number:	G380
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human LSHR. AA range:621-670
Concentration:	1 mg/ml
Gene Name:	LHCGR
Protein Name:	Lutropin-choriogonadotropic hormone receptor
Other Name:	LHCGR; LCGR; LGR2; LHRHR; Lutropin-choriogonadotropic hormone receptor; LH/CG-R; Luteinizing hormone receptor; LHR; LSH-R
Human Swiss Prot No.:	P22888
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	LHR Polyclonal Antibody detects endogenous levels of LHR protein.
Storage and Stability:	-20°C/1 year