

Antibodies Product Datasheet

Product Name:	Mas1-like Antibody
Catalog Number:	G383
Reactivity:	Human, Rat, Mouse
Applications:	IF/ICC, ELISA
Dilutions:	Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAS1L. AA range:31-80
Concentration:	1 mg/ml
Gene Name:	MAS1L
Protein Name:	Mas-related G-protein coupled receptor MRG
Other Name:	MAS1L; MRG; Mas-related G-protein coupled receptor MRG; MAS-R; MAS1-like
Human Swiss Prot No.:	P35410
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Mas1-like Polyclonal Antibody detects endogenous levels of Mas1-like protein.
Storage and Stability:	-20°C/1 year