

Antibodies Product Datasheet

Product Name:	mPR α Antibody
Catalog Number:	G388
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MPRA. AA range:297-346
Concentration:	1 mg/ml
Gene Name:	PAQR7
Protein Name:	Membrane progesterin receptor alpha
Other Name:	PAQR7; MRPA; Membrane progesterin receptor alpha; mPR alpha; Progesterin and adipoQ receptor family member 7; Progesterin and adipoQ receptor family member VII
Human Swiss Prot No.:	Q86WK9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	mPR α Polyclonal Antibody detects endogenous levels of mPR α protein.
Storage and Stability:	-20°C/1 year