

Antibodies Product Datasheet

Product Name:	MRGE Antibody
Catalog Number:	G391
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/40000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRGRE. AA range:171-220
Concentration:	1 mg/ml
Gene Name:	MRGPRE
Protein Name:	Mas-related G-protein coupled receptor member E
Other Name:	MRGPRE; GPR167; MRGE; Mas-related G-protein coupled receptor member E; G-protein coupled receptor 167
Human Swiss Prot No.:	Q86SM8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRGE Polyclonal Antibody detects endogenous levels of MRGE protein.
Storage and Stability:	-20°C/1 year