

Antibodies Product Datasheet

Product Name:	MRGF Antibody
Catalog Number:	G392
Reactivity:	Human, Mouse, Rat
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRGRF. AA range:271-320
Concentration:	1 mg/ml
Gene Name:	MRGPRF
Protein Name:	Mas-related G-protein coupled receptor member F
Other Name:	MRGPRF; GPR140; GPR168; MRGF; PSEC0142; Mas-related G-protein coupled receptor member F; Mas-related gene F protein; G-protein coupled receptor 140; G-protein coupled receptor 168
Human Swiss Prot No.:	Q96AM1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRGF Polyclonal Antibody detects endogenous levels of MRGF protein.
Storage and Stability:	-20°C/1 year