

Antibodies Product Datasheet

Product Name:	GPR38 Antibody
Catalog Number:	G397
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from the C-terminal region of human GPR38.
Concentration:	1 mg/ml
Gene Name:	MLNR
Protein Name:	Motilin receptor
Other Name:	MLNR; GPR38; MTLR; MTLR1; Motilin receptor; G-protein coupled receptor 38
Human Swiss Prot No.:	O43193
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR38 Polyclonal Antibody detects endogenous levels of GPR38 protein.
Storage and Stability:	-20°C/1 year