

Antibodies Product Datasheet

Product Name:	Olfactory receptor 52E6 Antibody
Catalog Number:	G462
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human OR52E6. AA range:264-313
Concentration:	1 mg/ml
Gene Name:	OR52E6
Protein Name:	Olfactory receptor 52E6
Other Name:	OR52E6; Olfactory receptor 52E6; Olfactory receptor OR11-58
Human Swiss Prot No.:	Q96RD3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Olfactory receptor 52E6 Polyclonal Antibody detects endogenous levels of Olfactory receptor 52E6 protein.
Storage and Stability:	-20°C/1 year