

Antibodies Product Datasheet

Product Name:	GPR172B Antibody
Catalog Number:	G723
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PEVR2. AA range:235-284
Concentration:	1 mg/ml
Gene Name:	SLC52A1
Protein Name:	Solute carrier family 52 riboflavin transporter member 1
Other Name:	SLC52A1; GPR172B; PAR2; RFT1; Solute carrier family 52; riboflavin transporter, member 1; Porcine endogenous retrovirus A receptor 2; PERV-A receptor 2; Protein GPR172B; Riboflavin transporter 1; hRFT1 Q9NWF4
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GPR172B Polyclonal Antibody detects endogenous levels of GPR172B protein.
Storage and Stability:	-20°C/1 year