

Antibodies Product Datasheet

Product Name:	T2R14 Antibody
Catalog Number:	G753
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human TAS2R14. AA range:229-278
Concentration:	1 mg/ml
Gene Name:	TAS2R14
Protein Name:	Taste receptor type 2 member 14
Other Name:	TAS2R14; Taste receptor type 2 member 14; T2R14; Taste receptor family B member 1; TRB1
Human Swiss Prot No.:	Q9NYV8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	T2R14 Polyclonal Antibody detects endogenous levels of T2R14 protein.
Storage and Stability:	-20°C/1 year