

Antibodies Product Datasheet

Product Name:	AVP Receptor V2 Antibody
Catalog Number:	G788
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human AVPR2. AA range:72-121
Concentration:	1 mg/ml
Gene Name:	AVPR2
Protein Name:	Vasopressin V2 receptor
Other Name:	AVPR2; ADHR; DIR; DIR3; V2R; Vasopressin V2 receptor; V2R; AVPR V2; Antidiuretic hormone receptor; Renal-type arginine vasopressin receptor
Human Swiss Prot No.:	P30518
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AVP Receptor V2 Polyclonal Antibody detects endogenous levels of AVP Receptor V2 protein.
Storage and Stability:	-20°C/1 year