

Antibodies Product Datasheet

Product Name:	MRGX3 Antibody
Catalog Number:	G816
Reactivity:	Human, Rat, Mouse
Applications:	WB, IF/ICC, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MRGX3. AA range:99-148
Concentration:	1 mg/ml
Gene Name:	MRGPRX3
Protein Name:	Mas-related G-protein coupled receptor member X3
Other Name:	MRGPRX3; MRGX3; SNSR1; SNSR2; Mas-related G-protein coupled receptor member X3; Sensory neuron-specific G-protein coupled receptor 1/2
Human Swiss Prot No.:	Q96LB0
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MRGX3 Polyclonal Antibody detects endogenous levels of MRGX3 protein.
Storage and Stability:	-20°C/1 year