

Antibodies Product Datasheet

Product Name:	MAGE-A5 Antibody
Catalog Number:	G875
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human MAGEA5. AA range:68-117
Concentration:	1 mg/ml
Gene Name:	MAGEA5
Protein Name:	Melanoma-associated antigen 5
Other Name:	MAGEA5; MAGE5; Melanoma-associated antigen 5; Cancer/testis antigen 1.5; CT1.5; MAGE-5 antigen
Human Swiss Prot No.:	P43359
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	MAGE-A5 Polyclonal Antibody detects endogenous levels of MAGE-A5 protein.
Storage and Stability:	-20°C/1 year