

## Antibodies Product Datasheet

|                               |                                                                                                                                                             |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | Cleaved-Caspase-9 (D353) Antibody                                                                                                                           |
| <b>Catalog Number:</b>        | L0110                                                                                                                                                       |
| <b>Reactivity:</b>            | Mouse,Rat                                                                                                                                                   |
| <b>Applications:</b>          | IF/ICC, WB, IHC-p, ELISA                                                                                                                                    |
| <b>Dilutions:</b>             | IF: 1:50-200 Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.                       |
| <b>Modification:</b>          | Cleaved Specific                                                                                                                                            |
| <b>Source:</b>                | Polyclonal, Rabbit                                                                                                                                          |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                       |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from human Caspase 9. AA range:323-372                                                       |
| <b>Concentration:</b>         |                                                                                                                                                             |
| <b>Gene Name:</b>             | CASP9                                                                                                                                                       |
| <b>Protein Name:</b>          | Caspase9                                                                                                                                                    |
| <b>Other Name:</b>            | CASP9; MCH6; Caspase-9; CASP-9; Apoptotic protease Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6; ICE-LAP6           |
| <b>Human Swiss Prot No.:</b>  |                                                                                                                                                             |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                     |
| <b>Specificity:</b>           | Cleaved-Caspase-9 (D353) Polyclonal Antibody detects endogenous levels of fragment of activated Caspase-9 protein resulting from cleavage adjacent to D353. |
| <b>Storage and Stability:</b> | Store at -20° C/1 year                                                                                                                                      |