

Antibodies Product Datasheet

Product Name:	ICAD Antibody
Catalog Number:	L0113
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human DFFA. AA range:175-224
Concentration:	1 mg/ml
Gene Name:	DFFA
Protein Name:	DNA fragmentation factor subunit alpha
Other Name:	DFFA; DFF1; DFF45; H13; DNA fragmentation factor subunit alpha; DNA fragmentation factor 45 kDa subunit; DFF-45; Inhibitor of CAD; ICAD
Human Swiss Prot No.:	O00273
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ICAD Polyclonal Antibody detects endogenous levels of ICAD protein.
Storage and Stability:	-20°C/1 year