

Antibodies Product Datasheet

Product Name:	Cleaved-Lamin A (D230) Antibody
Catalog Number:	L0116
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. Immunohistochemistry: 1/100 - 1/300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human Lamin A. AA range:181-230
Concentration:	1 mg/ml
Gene Name:	LMNA
Protein Name:	Prelamin-A/C
Other Name:	LMNA; LMN1; Prelamin-A/C
Human Swiss Prot No.:	P02545
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Lamin A (D230) Polyclonal Antibody detects endogenous levels of fragment of activated Lamin A protein resulting from cleavage adjacent to D230.
Storage and Stability:	-20°C/1 year