

Antibodies Product Datasheet

Product Name:	Acrosin Antibody
Catalog Number:	L0125
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300. Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Acrosin . at AA range: 40-120
Concentration:	1 mg/ml
Gene Name:	ACR
Protein Name:	Acrosin
Other Name:	ACR; ACRS; Acrosin
Human Swiss Prot No.:	P10323
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Acrosin Polyclonal Antibody detects endogenous levels of Acrosin protein.
Storage and Stability:	-20°C/1 year