

Antibodies Product Datasheet

Product Name:	ADAM10 Antibody
Catalog Number:	L0127
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ADAM10 . at AA range: 170-250
Concentration:	1 mg/ml
Gene Name:	ADAM10
Protein Name:	Disintegrin and metalloproteinase domain-containing protein 10
Other Name:	ADAM10; KUZ; MADM; Disintegrin and metalloproteinase domain-containing protein 10; ADAM 10; CDw156; Kuzbanian protein homolog; Mammalian disintegrin-metalloprotease; CD antigen CD156c
Human Swiss Prot No.:	O14672
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ADAM10 Polyclonal Antibody detects endogenous levels of ADAM10 protein.
Storage and Stability:	-20°C/1 year