

Antibodies Product Datasheet

Product Name:	Cleaved-Cathepsin L2 (L114) Antibody
Catalog Number:	L0185
Reactivity:	Human, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human CATL2. AA range:95-144
Concentration:	1 mg/ml
Gene Name:	CTSL2
Protein Name:	Cathepsin L2
Other Name:	CTSL2; CATL2; CTSU; CTSV; Cathepsin L2; Cathepsin U; Cathepsin V
Human Swiss Prot No.:	O60911
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Cathepsin L2 (L114) Polyclonal Antibody detects endogenous levels of fragment of activated Cathepsin L2 protein resulting from cleavage adjacent to L114.
Storage and Stability:	-20°C/1 year