

Antibodies Product Datasheet

Product Name:	HABP2 Antibody
Catalog Number:	L0272
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from HABP2 . at AA range: 270-350
Concentration:	1 mg/ml
Gene Name:	HABP2
Protein Name:	Hyaluronan-binding protein 2
Other Name:	HABP2; HGFAL; PHBP; Hyaluronan-binding protein 2; Factor VII-activating protease; Factor seven-activating protease; FSAP; Hepatocyte growth factor activator-like protein; Plasma hyaluronan-binding protein
Human Swiss Prot No.:	Q14520
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HABP2 Polyclonal Antibody detects endogenous levels of HABP2 protein.
Storage and Stability:	-20°C/1 year