

Antibodies Product Datasheet

Product Name:	Cleaved-Notch 2 (A1734) Antibody
Catalog Number:	L0352
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NOTCH2. AA range:1715-1764
Concentration:	1 mg/ml
Gene Name:	NOTCH2
Protein Name:	Neurogenic locus notch homolog protein 2
Other Name:	NOTCH2; Neurogenic locus notch homolog protein 2; Notch 2; hN2
Human Swiss Prot No.:	Q04721
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Notch 2 (A1734) Polyclonal Antibody detects endogenous levels of fragment of activated Notch 2 protein resulting from cleavage adjacent to A1734.
Storage and Stability:	-20°C/1 year