

Antibodies Product Datasheet

Product Name:	Cleaved-Notch 4 (V1432) Antibody
Catalog Number:	L0363
Reactivity:	Human, Mouse, Monkey
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/5000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human NOTCH4. AA range:1401-1450
Concentration:	1 mg/ml
Gene Name:	NOTCH4
Protein Name:	Neurogenic locus notch homolog protein 4
Other Name:	NOTCH4; INT3; Neurogenic locus notch homolog protein 4; Notch 4; hNotch4
Human Swiss Prot No.:	Q99466
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Notch 4 (V1432) Polyclonal Antibody detects endogenous levels of fragment of activated Notch 4 protein resulting from cleavage adjacent to V1432.
Storage and Stability:	-20°C/1 year