

Antibodies Product Datasheet

Product Name:	Cleaved-Plasminogen HC A short form (V98) Antibody
Catalog Number:	L0370
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Cleaved Specific
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from human PLMN. AA range:79-128
Concentration:	1 mg/ml
Gene Name:	PLG
Protein Name:	Plasminogen
Other Name:	PLG; Plasminogen
Human Swiss Prot No.:	P00747
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cleaved-Plasminogen HC A short form (V98) Polyclonal Antibody detects endogenous levels of fragment of activated Plasminogen HC A short form protein resulting from cleavage adjacent to V98.
Storage and Stability:	-20°C/1 year